

Air Express™ APM

Alternating Pressure Therapy



Model Number: P1000D



Air Express™ APM

Dynamic Alternating Pressure Mattress

The Air Express APM control unit is adjustable to the patients weight, extremely quiet and lightweight.



The Control Unit

- Dynamic alternating pressure every five minutes for optimum therapy.
- Simple dial adjustment for customized therapy.
- Provides pressure management throughout entire alternating cycle.
- Quick hose release for rapid CPR deflation.
- Compact and lightweight.

The Therapeutic Support Surface

- Dynamic Alternating Pressure every five minutes.
- All air cushions are constructed of high quality and very pliable coated nylon for low shear and even lower interface pressures.
- Convoluted foam base provides additional safety to your patient.
- The PLEXUS® Quilted Top Sheet is highly vapor permeable and oversized providing a low friction, low shear barrier against complications due to moisture. Machine washable.
- Modular design facilitates interchangeability, cleaning and setup.
- Effective therapy for all stages of pressure ulcers.

FootStop/ZIP™ Cushion (Optional)

- Interchangeable Foot Stop/ZIP™ cushion permits tailoring the mattress length to fit the patient while reducing pressure in the critical heel area.
- ZIP cushion can be used independently anywhere along the mattress.

The optional Foot Stop/ZIP™ cushion uniquely maintains proper patient positioning while *reducing* interface pressure on the critical heel area.



The therapeutic support surface is lightweight and made of high quality materials. It conveniently rolls up for easy carrying with attached straps. Separate carrying case also available.



Specifications: P1000D Mattress Replacement

| | |
|----------------------|---|
| Dimensions: | Control Unit: 5" x 10" x 12" / 12.7cm x 25.4cm x 30.5cm (5lbs. / 2.3kg) Mattress: 10" x 34" x 80" / 25.4cm x 86.3cm x 203.2cm (15lbs. / 6.8kg) |
| Electrical Rating: | 120V, 60Hz, 1A / 230V, 50Hz, .5A |
| Agency Approvals: | UL 2601/CAN/CSA-C22.2 No. 601-1-M90, TUV, EN 60601-1 |
| AP Cycle Time: | 5 Minutes |
| Warranty: | *Control Unit: 2 year Ltd. Mattress: 1 year Ltd. |
| Weight Limit: | 350 lbs. |
| Medicare HCPCS Code: | E0277 |

*P1000I Economy Control Unit available - 1 year warranty

*Available in 220 Volt model

Available in Canada through:
Waterloo Bedding Company Ltd.
825 Trillium Drive
Kitchener, Ontario N2R 1J9
Phone: 1-800-203-4293 • 519-742-4447
Fax: 519-742-3699



Plexus Medical • A Gaymar Company
10 Centre Drive, Orchard Park, NY 14127
Toll Free 1.800.828.7341
Fax Free 1.800.993.7890
International 716.662.8636
International Fax 716.662.0730
www.gaymar.com

GAYMAR®

MAW5111 • 09/03